

Naive CD4+T

Catalog No.: RK50042

Basic Information

Observed MW
Sanger Sequencing

Calculated MW

Category
Primary antibody

Applications

Cross-Reactivity

Background

Recommended Dilutions

Immunogen Information

Gene ID	Swiss Prot
Immunogen	
Synonyms	

Contact

	400-999-6126
	cn.market@abclonal.com.cn
	www.abclonal.com.cn

Product Information

Source	Isotype	Purification
Mouse		
Storage		
Store at 2-8°C.		