

ZEB1 Knockout HeLa Cell Line, Homozygous

Catalog No.: RM02165

Basic Information

Catalog No.

RM02165

Category

Cell Line

Parental Cell line

HeLa

Genotype

Knockout

Gene Information

Gene Symbol

ZEB1

Species

Human

Gene ID

6935

Swiss Prot

P37275

Synonyms

AREB6; BZP; DELTAEF1; FECD6; NIL2A;
PPCD3; TCF8; ZFHEP; ZFH1A

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Background

This gene encodes a zinc finger transcription factor. The encoded protein likely plays a role in transcriptional repression of interleukin 2. Mutations in this gene have been associated with posterior polymorphous corneal dystrophy-3 and late-onset Fuchs endothelial corneal dystrophy. Alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Mar 2010]

Product Information

Description

ZEB1 Knockout HeLa Cell Line is engineered from HeLa cell line with Gene-Editing Technology.

Allele-1:64bp deletion in exon6

Allele-2:64bp deletion in exon6

Mammalian cells such as human, rat and mouse cells are normally diploid with two alleles. Homozygote: both alleles were knocked out, mRNA has no signal, no expression of proteins. Heterozygote: only one allele was knocked out, the mRNA transcript levels was decreased compared to wild type, and the protein expression levels was also lower than that of the wild type.

Packaging

1 vial parental cell line and 1 vial knockout cell line

Shipping Conditions

Dry ice

Amount

1~5x10⁶ cells/vial

Storage

Stored in liquid nitrogen for a long time less than -130°C. Minimizing freeze-thaw cycles.

Protocol

Upon arrival, it should be maintained in DMEM medium with 10%(v/v) fetal bovine serum and 100U penicillin-streptomycin, at 37°C with 5% CO₂ condition.

1. Thaw the vial in 37°C water bath, and shake it to melt as soon as possible.
2. Transfer the cell suspension to a 15mL conical tube with pre-warmed 5mL complete medium and centrifuge 1000rpm for approximately 5 minutes at room temperature.
3. Remove and discard the supernatant.
4. Resuspend the cell pellet with 1mL pre-warmed complete medium and seed in 10cm dish.
5. Add 8-10mL of complete medium.
6. Incubate the culture at 37°C incubator with 5% CO₂.
7. A subcultivation ratio of 1:2-1:4 is recommended.

Sequencing data

WT AATGGGCGACCAAG*****CGACCACAGATACG
Mut AATGGGCGACCAAG***Deletion***CGACCACAGATACG
Allele-1: 64bp deletion in exon6
WT AATGGGCGACCAAG*****CGACCACAGATACG
Mut AATGGGCGACCAAG***Deletion***CGACCACAGATACG
Allele-2: 64bp deletion in exon6

Genome sequence analysis of PCR products from parental (WT) and ZEB1 knockout (KO) HeLa cells, using sanger sequencing.