

EEF2K Knockout NIH/3T3 Cell Lysate, Homozygous

Catalog No.: RM02331

Basic Information

Catalog No.

RM02331

Category

Cell Lysate

Parental Cell line

NIH/3T3

Genotype

Knockout

Gene Information

Gene Symbol

EEF2K

Species

Mouse

Gene ID

13631

Synonyms

C86191; eEF-2K

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Background

Product Information

Description

EEF2K Knockout NIH/3T3 Cell Line is engineered from NIH/3T3 cell line with Gene-Editing technology.

Allele-1:1bp deletion in exon1

Allele-2:1bp deletion in exon1

Mammalian cells such as human, rat and mouse cells are normally diploid with two alleles. Homozygote: both alleles were knocked out, mRNA has no signal, no expression of proteins. Heterozygote: only one allele was knocked out, the mRNA transcript levels was decreased compared to wild type, and the protein expression levels was also lower than that of the wild type.

Packaging

1 vial parental cell Lysate and 1 vial knockout cell Lysate

Shipping Conditions

4°C

Amount

50µL, 2µg/µL.

Storage

Lysate is stable for 12 months when stored at -20°C. Minimizing freeze-thaw cycles.

Protocol

To be used as WB control. Lysate is supplied in 1× SDS sample buffer (2% SDS, 60 mM Tris-HCl pH 6.8, 10% Glycerol, 0.02% Bromophenol blue, 60 mM beta-mercaptoethanol). Lysate should be boiled for 3 - 5 minutes before loading onto gel.

Sequencing data

WT CAGAGTGC GGCTCCACAGGGTCACCAGCC
Mut CAGAGTGC GGCTCC -CAGGGTCACCAGCC
Allele-1: 1bp deletion in exon1

WT CAGAGTGC GGCTCCACAGGGTCACCAGCC
Mut CAGAGTGC GGCTCC -CAGGGTCACCAGCC
Allele-2: 1bp deletion in exon1

Genome sequence analysis of PCR products from parental (WT) and EEF2K knockout (KO) NIH3T3 cells, using sanger sequencing.