

FAM50A Knockout NIH/3T3 Cell Lysate, Homozygous

Catalog No.: RM02338

Basic Information

Catalog No.

RM02338

Category

Cell Lysate

Parental Cell line

NIH/3T3

Genotype

Knockout

Gene Information

Gene Symbol

FAM50A

Species

Mouse

Gene ID

108160

Synonyms

D0HXS9928E; XAP-5

Contact

<u>a</u>	400-	999-6126
\bowtie	cn.market@abclona	al.com.cn
•	www.abclona	al.com.cn

Background

Product Information

Description

FAM50A Knockout NIH/3T3 Cell Line is engineered from NIH/3T3 cell line with Gene-Editing technology.

Allele-1:exon2 was deleted

Allele-2:exon2 was deleted

Mammalian cells such as human, rat and mouse cells are normally diploid with two alleles. Homozygote: both alleles were knocked out, mRNA has no signal, no expression of proteins. Heterozygote: only one allele was knocked out, the mRNA transcript levels was decreased compared to wild type, and the protein expression levels was also lower than that of the wild type.

Packaging

 ${\bf 1}$ vial parental cell Lysate and ${\bf 1}$ vial knockout cell Lysate

 $\begin{array}{ll} \textbf{Shipping Conditions} & \textbf{Amount} \\ 4^{\circ} C & 50 \mu L, 2 \mu g/\mu L. \end{array}$

Storage

Lysate is stable for 12 months when stored at -20°C. Minimizing freeze-thaw cycles.

Protoco

To be used as WB control. Lysate is supplied in $1 \times$ SDS sample buffer (2% SDS, 60 mM Tris-HCl pH 6.8, 10% Glycerol, 0.02% Bromophenol blue, 60 mM beta-mercaptoethanol). Lysate should be boiled for 3 - 5 minutes before loading onto gel.

Sequencing data

WT GTTAAGGCA*TCGA***********TGTCACTCTTTTGC
Mut GTTA-GGCA*TCGA***Deletion***TGTCACTCTTTTGC
Allele-1: exon2 was deleted

WT GTTAAGGCA*TCGA*************TGTCACTCTTTTGC
Mut GTTA-GGCA*TCGA***Deletion***TGTCACTCTTTTGC
Allele-2: exon2 was deleted

Genome sequence analysis of PCR products from parental (WT) and FAM50A Knockout (KO) NIH3T3 cells, using sanger sequencing.