

SALL4 Knockdown HeLa Cell Line, Heterozygous

Catalog No.: RM02577

Basic Information

Catalog No.

RM02577

Category

Cell Line

Parental Cell line

HeLa

Genotype

Knockdown

Gene Information

Gene Symbol

SALL4

Species

Human

Gene ID

57167

Swiss Prot

Q9UJQ4

Synonyms

DRRS; HSAL4; ZNF797

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Background

This gene encodes a zinc finger transcription factor thought to play a role in the development of abducens motor neurons. Defects in this gene are a cause of Duane-radial ray syndrome (DRRS). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015]

Product Information

Description

SALL4 Knockdown HeLa Cell Line is engineered from HeLa cell line with Gene-Editing technology.

Allele-1:144bp deletion in exon2

Allele-2:143bp deletion in exon2

Mammalian cells such as human, rat and mouse cells are normally diploid with two alleles. Homozygote: both alleles were knocked out, mRNA has no signal, no expression of proteins. Heterozygote: only one allele was knocked out, the mRNA transcript levels was decreased compared to wild type, and the protein expression levels was also lower than that of the wild type.

Packaging

1 vial parental cell line and 1 vial knockout cell line

Shipping Conditions

Dry ice

Amount

1~5x10⁶ cells/vial

Storage

Stored in liquid nitrogen for a long time less than -130°C. Minimizing freeze-thaw cycles.

Protocol

Upon arrival, it should be maintained in DMEM medium with 10%(v/v) fetal bovine serum and 100U penicillin-streptomycin, at 37°C with 5% CO₂ condition.

1. Thaw the vial in 37°C water bath ,and shake it to melt as soon as possible.
2. Transfer the cell suspension to a 15mL conical tube with pre-warmed 5mL complete medium and centrifuge 1000rpm for approximately 5 minutes at room temperature.
3. Remove and discard the supernatant.
4. Resuspend the cell pellet with 1mL pre-warmed complete medium and seed in 10cm dish.
5. Add 8-10mL of complete medium.
6. Incubate the culture at 37°C incubator with 5% CO₂.
7. A subcultivation ratio of 1:2-1:4 is recommended.

Sequencing data

WT GGAGAAGCCGGATG*****GCGCGGATGCACTC
Mut GGAGAAGCCGGATG***Deletion***GCGCGGATGCACTC
Allele-1: 144bp deletion in exon2

WT GGAGAAGCCGGATG*****AGCGCGGATGCACT
Mut GGAGAAGCCGGATG***Deletion***AGCGCGGATGCACT
Allele-2: 143bp deletion in exon2

Genome sequence analysis of PCR products from parental (WT) and SALL4 Knockdown (KD) HeLa cells, using sanger sequencing.