

# Rabbit anti-Human IL-6 mAb (DET)

**Catalog No.: RM17602**

## Basic Information

**Catalog No.**  
RM17602

**Catagory**  
Elisa Antibody Kit

**Application**  
ELISA

## Product Information

**Ig Type**  
Rabbit IgG

**Purification**  
Affinity purification

**Endotoxin Level**  
< 1.0 EU/µg of the protein by LAL method.

**Storage**  
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2um filtered solution in PBS with 0.1%Braveds MB-1,PH 7.4.

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

This gene encodes a cytokine that functions in inflammation and the maturation of B cells. In addition, the encoded protein has been shown to be an endogenous pyrogen capable of inducing fever in people with autoimmune diseases or infections. The protein is primarily produced at sites of acute and chronic inflammation, where it is secreted into the serum and induces a transcriptional inflammatory response through interleukin 6 receptor, alpha. The functioning of this gene is implicated in a wide variety of inflammation-associated disease states, including susceptibility to diabetes mellitus and systemic juvenile rheumatoid arthritis. Elevated levels of the encoded protein have been found in virus infections, including COVID-19 (disease caused by SARS-CoV-2).

## Immunogen Information

**Immunogen**  
Recombinant IL-6 Protein

## Cross-Reactivity

No cross-reactivity in ELISA assay with recombinant hCntf, hG-CSF, hGM-CSF, hIL-1A, hIL-1B, hIL-2, hIL-3, hIL-4, hIL-6RA, hIL-6RB, hIL-7, hIL-8, hIL-11, hIL-12, hILF, hOSM, hTNFA, mL-2, mL-4, mL-5, mL-6, mL-11, mL-12


## Assay Applications

Human IL-6 Sanwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4ug/mL	Rabbit anti-Human IL-6(CAP)(Cat. No.RM17601)
ELISA Detection	0.033-0.13ug/mL	Rabbit anti-Human IL-6(DET)(Cat. No.RM17602)
Standard	1.56-100pg/mL	Recombinant Human IL-6 Protein

## Validation Data

---



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.