

Rabbit anti-Human NEFL mAb (DET)

Catalog No.: RM17724

Basic Information

Catalog No.
RM17724

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification


Endotoxin Level


Storage
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS with 0.1%Braveds MB-1,PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Neurofilaments usually contain three intermediate filament proteins: NEFL, NEFM, and NEFH which are involved in the maintenance of neuronal caliber. May additionally cooperate with the neuronal intermediate filament proteins PRPH and INA to form neuronal filamentous networks (By similarity).

Immunogen Information

Immunogen
Recombinant Human NEFL Protein


Cross-Reactivity

Assay Applications

Human NEFL Sanwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	0.5-2ug/mL	Rabbit anti-Human NEFL(CAP)(Cat. No.RM17718)
ELISA Detection	0.025-0.1ug/mL	Rabbit anti-Human NEFL(DET)(Cat. No.RM17719)
Standard	31.25-2000pg/mL	Recombinant Human NEFL Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.