

# Rabbit anti-Human LGALS3 mAb (DET)

**Catalog No.:** RM17862

## Basic Information

**Catalog No.**  
RM17862

**Catagory**  
Elisa Antibody Kit

**Application**  
ELISA

## Product Information

**Ig Type**  
Rabbit IgG


**Purification**  
Affinity purification


**Endotoxin Level**


**Storage**  
Store at -20°C. Avoid freeze / thaw cycles. Preservative[]0.05%ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2um filtered solution in PBS with 0.05%ProClin 300,PH 7.4.

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

This gene encodes a member of the galectin family of carbohydrate binding proteins. Members of this protein family have an affinity for beta-galactosides. The encoded protein is characterized by an N-terminal proline-rich tandem repeat domain and a single C-terminal carbohydrate recognition domain. This protein can self-associate through the N-terminal domain allowing it to bind to multivalent saccharide ligands. This protein localizes to the extracellular matrix, the cytoplasm and the nucleus. This protein plays a role in numerous cellular functions including apoptosis, innate immunity, cell adhesion and T-cell regulation. The protein exhibits antimicrobial activity against bacteria and fungi. Alternate splicing results in multiple transcript variants.

## Immunogen Information

**Immunogen**  
Recombinant Human LGALS3

## Cross-Reactivity

/


## Assay Applications

Human LGALS3 Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4ug/mL	Rabbit anti-Human LGALS3 (CAP)(Cat. No.RM17861)
ELISA Detection	0.027-0.111ug/mL	Rabbit anti-Human LGALS3 (DET)(Cat. No.RM17862)
Standard	0.15625-10ng/mL	Recombinant Human LGALS3 Protein

## Validation Data

---



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.