

# Rabbit anti-Human uPAR mAb(DET)

**Catalog No.:** RM18034

## Basic Information

**Catalog No.**  
RM18034

**Catagory**  
Elisa Antibody Kit

**Application**  
ELISA

## Product Information

**Ig Type**  
Rabbit IgG


**Purification**  
Affinity purification


**Endotoxin Level**

**Storage**  
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.  
Preservative 0.05% ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2µm filtered solution in PBS with 0.05% ProClin 300, pH 7.4 containing human uPAR Antibody.

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

Acts as a receptor for urokinase plasminogen activator. Plays a role in localizing and promoting plasmin formation. Mediates the proteolysis-independent signal transduction activation effects of U-PA. It is subject to negative-feedback regulation by U-PA which cleaves it into an inactive form.

## Immunogen Information

**Immunogen**  
Recombinant Human uPAR

## Cross-Reactivity

/


## Assay Applications

Human uPAR Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4µg/mL	Rabbit anti-Human uPAR (CAP)(Cat. No.RM18033)
ELISA Detection	0.05-0.2µg/mL	Rabbit anti-Human uPAR (DET)(Cat. No.RM18034)
Standard	62.5-4000pg/mL	Recombinant Human uPAR Protein

## Validation Data

---



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.