

Rabbit anti-Human Ghrelin mAb(CAP)

Catalog No.: RM18049

Basic Information

Catalog No.
RM18049

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.


Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS with 0.1% BSA, pH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Ghrelin is the ligand for growth hormone secretagogue receptor type 1 (GHSR) Induces the release of growth hormone from the pituitary Has an appetite-stimulating effect, induces adiposity and stimulates gastric acid secretion. Involved in growth regulation

Immunogen Information

Immunogen
Recombinant Human Ghrelin

Cross-Reactivity


/

Assay Applications

Human Ghrelin Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	2-8ug/mL	Rabbit anti-Human Ghrelin mAb(CAP)(Cat. No.RM18049)
ELISA Detection	0.008-0.032ug/mL	Rabbit anti-Human Ghrelin mAb(DET)(Cat. No.RM18050)
Standard	0.16-10ng/mL	Recombinant Human Ghrelin Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.