www.abclonal.com

ABclonal

Rabbit anti-Mouse/Rat IFN-gamma mAb

Catalog No.: RMK0001

Basic Information

Catalog No.

RMK0001

Catagory

Elisa Antibody Kit

Application

ELISA, Multiplex

Product Information

Ig Type

Rabbit IgG

Purification

Affinity purification

Endotoxin Level

/

Storage

Store at -20°C.

Avoid repeated freeze-thaw cycles.

Formulation

Supplied as a 0.2µm filtered solution in PBS with 0.05% Proclin300, PH 7.4.

Contact

•	order@abclonal.com
<u>a</u>	support@abclonal.com
•	www.abclonal.com

Background

Interferon-gamma (IFN-gamma), also known as type II or immune interferon, exerts a wide range of immunoregulatory activities and is considered to be the prototype proinflammatory cytokine . Mature rat IFN-gamma exists as a noncovalently linked homodimer of 20 - 25 kDa variably glycosylated subunits. It shares 86% amino acid sequence identity with mouse IFN-gamma and 37%?45% with bovine, canine, cotton rat, equine, feline, human, porcine, and rhesus IFN-gamma. IFN-gamma dimers bind to IFN-gamma RI (alpha subunits) which then interact with IFN-gamma RII (beta subunits) to form the functional receptor complex of two alpha and two beta subunits. Inclusion of IFN-gamma RII increases the binding affinity for ligand and the efficiency of signal transduction. IFN? gamma is produced by a variety of immune cells under inflammatory conditions, notably by T cells and NK cells. It plays a key role in host defense by promoting the development and activation of Th1 cells, chemoattraction and activation of monocytes and macrophages, upregulation of ...

Immunogen Information

Immunogen

Recombinant Rat IFN-gamma Protein.

Cross-Reactivity


Mouse, Rat

Assay Applications


Mouse/Rat IFN-gamma Sandwich Immunoassay

		Recommended Concentration	Sample
ELISA	Capture	1-4 μg/mL	Rabbit anti-Mouse/Rat IFN-gamma mAb (Cat. No. RMK0001)
	Detection	0.05-0.2 μg/mL	Rabbit anti-Mouse IFN-gamma mAb (Cat. No. RMK0582)
	Standard	15.63-1000 pg/mL	Recombinant Mouse IFN-gamma Protein (Cat. No. RP01070)

Multiplex	Capture	3-20 μg/mL	Rabbit anti-Mouse/Rat IFN-gamma mAb (Cat. No. RMK0001)
	Detection	0.017-2 μg/mL	Rabbit anti-Rat IFN-gamma mAb (Cat. No.RMK0002)
	Standard	1.37-1000 pg/mL	Recombinant Rat IFN-gamma Protein
Multiplex	Capture	0.1-5 μg/ml	Rabbit anti-Mouse/Rat IFN-gamma mAb (Cat. No. RMK0001)
	Detection	0.5-2 μg/ml	Rabbit anti-Rat IFN-gamma mAb (Cat. No.RMK0002)
	Standard	6.86-5000 pg/mL	Recombinant Rat IFN-gamma Protein



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.