

Rabbit anti-Human IL-18 mAb(DET)

Catalog No.: RMK0032

Basic Information

Catalog No.
RMK0032

Catagory
Elisa Antibody Kit

Application
ELISA, Multiplex

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
/

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

The protein encoded by this gene is a proinflammatory cytokine of the IL-1 family that is constitutively found as a precursor within the cytoplasm of a variety of cells including macrophages and keratinocytes. The inactive IL-18 precursor is processed to its active form by caspase-1, and is capable of stimulating interferon gamma production, and of regulating both T helper (Th) 1 and Th2 responses. This cytokine has been implicated in the injury of different organs, and in potentially fatal conditions characterized by a cytokine storm. In humans, IL-18 gene is located on chromosome 11. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Immunogen Information

Immunogen
Recombinant Human IL-18 Protein

Cross-Reactivity


Assay Applications

Human IL-18 Sandwich Immunoassay

		Recommended Concentration	Sample
ELISA	Capture	1-4 µg/mL	Rabbit anti-Human IL-18 mAb (Cat. No. RMK0031)
	Detection	0.08-0.33 µg/mL	Rabbit anti-Human IL-18 mAb (Cat. No. RMK0032)
	Standard	15.6-1000 pg/mL	Recombinant Human IL-18 Protein

Multiplex	Capture	3-20 µg/ml	Rabbit anti-Human IL-18 mAb (Cat. No. RMK0031)
	Detection	0.017-2 µg/ml	Rabbit anti-Human IL-18 mAb (Cat. No. RMK0032)
	Standard	6.86-5000 pg/mL	Recombinant Human IL-18 Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.