

Rabbit anti-Human IL-12 p70 mAb

Catalog No.: RMK0070

Basic Information

Catalog No.
RMK0070

Category
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level


Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS with 0.05%ProClin 300,PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Immunogen Information

Immunogen
Recombinant Human IL-12P70

Cross-Reactivity


Assay Applications

Human IL-12P70 Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4ug/mL	Rabbit anti-Human IL-12P70 (CAP)(Cat. No.RMK0070)
ELISA Detection	0.05-0.2ug/mL	Rabbit anti-Human IL-12P70 (DET)(Cat. No.RM17740)
Standard	31.25-2000pg/mL	Recombinant Human IL-12P70 Protein


Validation Data

Human IL-12P70 ELISA Kit



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.

Stability



Placed at 37°C for 3 days, the stability of the standard curve all conform to CV <10%.