

Rabbit anti-Human S100A7 mAb(DET)

Catalog No.: RMK0080

Basic Information

Catalog No.
RMK0080

Category
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.
Preservative □ 0.05% ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2µm filtered solution in PBS with 0.05% ProClin 300, PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Immunogen Information

Immunogen
Recombinant Human S100A7

Cross-Reactivity

Assay Applications

Human S100A7 Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	2-8µg/mL	Rabbit anti-Human S100A7 (CAP)(Cat. No.RMK0079)
ELISA Detection	0.1-0.4µg/mL	Rabbit anti-Human S100A7 (DET)(Cat. No.RMK0080)
Standard	93.75-6000pg/mL	Recombinant Human S100A7 Protein