

# Rabbit anti-Human Serum Albumin mAb(DET)

**Catalog No.:** RMK0082

## Basic Information

**Catalog No.**  
RMK0082

**Catagory**  
Elisa Antibody Kit

**Application**  
ELISA

## Product Information

**Ig Type**  
Rabbit IgG


**Purification**  
Affinity purification


**Endotoxin Level**

**Storage**  
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.  
Preservative □ 0.05% ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2µm filtered solution in PBS with 0.05% ProClin 300, PH 7.4.

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

## Immunogen Information

**Immunogen**  
Native Human Serum albumin

## Cross-Reactivity

## Assay Applications

Human Serum Albumin Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	0.5-2µg/mL	Rabbit anti-Human Serum Albumin mAb(CAP)(Cat. No.RMK0081 )
ELISA Detection	0.016-0.066µg/mL	Rabbit anti-Human Serum Albumin mAb(DET)(Cat. No.RMK0082 )
Standard	1.56-100ng/mL	Recombinant Human Serum Albumin Protein