

Rabbit anti-Human GM-CSF mAb(CAP)

Catalog No.: RMK0105

Basic Information

Catalog No.
RMK0105

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
N/A

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of granulocytes and macrophages. The active form of the protein is found extracellularly as a homodimer. This gene has been localized to a cluster of related genes at chromosome region 5q31, which is known to be associated with interstitial deletions in the 5q- syndrome and acute myelogenous leukemia. Other genes in the cluster include those encoding interleukins 4, 5, and 13. This gene plays a role in promoting tissue inflammation. Elevated levels of cytokines, including the one produced by this gene, have been detected in SARS-CoV-2 infected patients that develop acute respiratory distress syndrome. Mice deficient in this gene or its receptor develop pulmonary alveolar proteinosis.

Immunogen Information

Immunogen
Recombinant Human GM-CSF Protein

Cross-Reactivity


No cross-reactivity in ELISA with recombinant H.G-CSF,H.IL-3,H.IL-4,M.IL-4,H.IL-5,M.IL-5,H.M-CSF

Assay Applications

Human GM-CSF Sandwich Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4 µg/mL	Rabbit anti-Human GM-CSF mAb (Cat. No. RMK0105)
ELISA Detection	0.033-0.132 µg/mL	Rabbit anti-Human GM-CSF mAb (Cat. No. RMK0107)
Standard	7.81-500 pg/mL	Recombinant Human GM-CSF Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.