

# Rabbit anti- Human TSLP mAb (CAP)

**Catalog No.: RMK0113**

## **Basic Information**

---

**Catalog No.**  
RMK0113

**Category**  
Elisa Antibody Kit

**Application**  
ELISA

## **Product Information**

---

**Ig Type**  
Rabbit IgG

**Purification**  
Affinity purification

**Endotoxin Level**


**Storage**  
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.


**Avoid repeated freeze-thaw cycles.**


**Formulation**  
Supplied as a 0.2um filtered solution in PBS with 0.1% BSA, pH 7.8.

## **Contact**

---

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

---

## **Background**

---

## **Immunogen Information**

---

**Immunogen**  
Recombinant Human TSLP

## **Cross-Reactivity**

---

## **Assay Applications**

---

Human TSLP Sandwich ELISA Immunoassay