

Rabbit anti-Mouse C-X-C motif chemokine 13 (CXCL13) mAb(DET)

Catalog No.: RMK0173

Basic Information

Catalog No.
RMK0173

Category
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.
Preservative 0.05% ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2µm filtered solution in PBS with 0.05% ProClin 300, PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Mouse CXCL13 (C-X-C Motif Chemokine Ligand 13), also known as B-Lymphocyte Chemoattractant (BLC), B-Cell Attracting Chemokine 1 (BCA-1), or Small Inducible Protein A14 (SIGLEC14), is a small secreted protein that belongs to the CXC subfamily of chemokines. Chemokines are a group of cytokines that play crucial roles in directing the migration of immune cells during inflammation and homeostasis. In mammals, CXCL13 is primarily expressed by follicular dendritic cells within germinal centers of secondary lymphoid organs such as spleen and lymph nodes. It functions as a potent chemoattractant for B cells expressing its receptor, CXCR5. The CXCL13/CXCR5 axis is particularly important for the organization of lymphoid tissue and the formation of B-cell follicles, which are specialized compartments where B cells undergo proliferation, differentiation, and affinity maturation in response to antigens. Research involving Mouse CXCL13 has been instrumental in understanding various aspects of humoral immunity and autoimmune diseases. Elevated levels...

Immunogen Information

Immunogen
Recombinant Mouse C-X-C motif chemokine 13 (CXCL13)

Cross-Reactivity

/

Assay Applications

Mouse C-X-C motif chemokine 13 (CXCL13) Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4µg/mL	Rabbit anti-Mouse C-X-C motif chemokine 13 (CXCL13) (CAP)(Cat. No.RMK0172)
ELISA Detection	0.05-0.2ng/mL	Rabbit anti-Mouse C-X-C motif chemokine 13 (CXCL13) (DET)(Cat. No.RMK0173)
Standard	31.25-2000pg/mL	Recombinant Mouse C-X-C motif chemokine 13 (CXCL13) Protein