

Rabbit anti-Human Complement Factor D mAb(CAP)

Catalog No.: RMK0176

Basic Information

Catalog No.
RMK0176

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
/

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This gene encodes a member of the S1, or chymotrypsin, family of serine peptidases. This protease catalyzes the cleavage of factor B, the rate-limiting step of the alternative pathway of complement activation. This protein also functions as an adipokine, a cell signaling protein secreted by adipocytes, which regulates insulin secretion in mice. Mutations in this gene underlie complement factor D deficiency, which is associated with recurrent bacterial meningitis infections in human patients. Alternative splicing of this gene results in multiple transcript variants. At least one of these variants encodes a preproprotein that is proteolytically processed to generate the mature protease.

Immunogen Information

Immunogen
Recombinant Human Complement Factor D

Cross-Reactivity

Human

Assay Applications

Human Complement Factor D Sandwich Immunoassay

	Recommended Concentration	Sample
ELISA Capture	0.5-2 µg/mL	Rabbit anti-Human Complement Factor D mAb (Cat. No. RMK0176)
ELISA Detection	0.033-0.133 µg/mL	Rabbit anti-Human Complement Factor D mAb (Cat. No. RMK0177)
Standard	7.81-500 pg/mL	Recombinant Human Complement Factor D Protein (Cat. No. RP01367)

Validation Data

