

Rabbit anti-Human Cadherin-3 (CDH3) mAb(CAP)

Catalog No.: RMK0185

Basic Information

Catalog No.
RMK0185

Category
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.


Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS with 0.1% BSA, pH 7.8.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Immunogen Information

Immunogen
Recombinant Human CDH3

Cross-Reactivity

/

Assay Applications

Human CDH3 Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4ug/mL	Rabbit anti-Human CDH3 mAb(CAP)(Cat. No.RMK0185)
ELISA Detection	0.0125-0.5ug/mL	Rabbit anti-Human CDH3 mAb(DET)(Cat. No.RMK0186)
Standard	62.5-4000pg/mL	Recombinant Human CDH3 Protein