

# Rabbit anti-Human Interleukin 6 receptor (IL-6R) mAb(DET)

**Catalog No.:** RMK0196

## Basic Information

**Catalog No.**  
RMK0196

**Catagory**  
Elisa Antibody Kit

**Application**  
ELISA

## Product Information

**Ig Type**  
Rabbit IgG


**Purification**  
Affinity purification


**Endotoxin Level**

**Storage**  
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2um filtered solution in PBS with 0.1%Braveds MB-1,PH 7.8.

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

This gene encodes a subunit of the interleukin 6 (IL6) receptor complex. Interleukin 6 is a potent pleiotropic cytokine that regulates cell growth and differentiation and plays an important role in the immune response. The IL6 receptor is a protein complex consisting of this protein and interleukin 6 signal transducer (IL6ST/GP130/IL6-beta), a receptor subunit also shared by many other cytokines. Dysregulated production of IL6 and this receptor are implicated in the pathogenesis of many diseases, such as multiple myeloma, autoimmune diseases and prostate cancer. Alternatively spliced transcript variants encoding distinct isoforms have been identified in this gene. A pseudogene of this gene is found on chromosome 9.

## Immunogen Information

**Immunogen**  
Recombinant Human IL-6R

## Cross-Reactivity

/

## Assay Applications

Human Interleukin 6 receptor (IL-6R) Sanwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4ug/mL	Rabbit anti-Human Interleukin 6 receptor (IL-6R) (CAP)(Cat. No.RMK0195)
ELISA Detection	0.05-0.2ng/mL	Rabbit anti-Human Interleukin 6 receptor (IL-6R) (DET)(Cat. No.RMK0196)
Standard	62.5-4000pg/mL	Recombinant Human Interleukin 6 receptor (IL-6R) Protein

## Validation Data

---

