

Rabbit anti-Mouse IL-13 mAb (CAP)

Catalog No.: RMK0202

Basic Information

Catalog No.
RMK0202

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification


Endotoxin Level


Storage
This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.


Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human IL-1 β Antibody.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

IL-13 (Interleukin-13) plays a key role in the pathogenesis of allergy, cancer, and tissue fibrosis. It is secreted by several helper T cell subsets, NK cells, mast cells, eosinophils, basophils, and visceral smooth muscle cells

Immunogen Information

Immunogen
Recombinant Mouse IL-13 Protein

Cross-Reactivity

/

Assay Applications

Mouse IL-13 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Mouse IL-13 mAb(Cat. No.RMK0202)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Mouse IL-13 mAb(Cat. No.RMK0203)
Standard	6.86-5000pg/mL	Recombinant Mouse IL-13 Protein