

Rabbit anti-Mouse VCAM-1/CD106 mAb(CAP)

Catalog No.: RMK0268

Basic Information

Catalog No.
RMK0268

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
/

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Predicted to enable integrin binding activity and primary amine oxidase activity. Acts upstream of or within several processes, including cellular response to glucose stimulus; chorio-allantoic fusion; and heterophilic cell-cell adhesion via plasma membrane cell adhesion molecules. Located in apical plasma membrane. Is expressed in several structures, including brain; embryo mesenchyme; extraembryonic component; heart; and hemolymphoid system. Human ortholog(s) of this gene implicated in hypertension. Orthologous to human VCAM1 (vascular cell adhesion molecule 1).

Immunogen Information

Immunogen
Recombinant mouse VCAM-1/CD106

Cross-Reactivity


/

Assay Applications

Mouse VCAM-1/CD106 Sandwich Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4 µg/mL	Rabbit anti-Mouse VCAM-1/CD106 mAb (Cat. No. RMK0268)
ELISA Detection	0.025-0.1 µg/mL	Rabbit anti-Mouse VCAM-1/CD106 mAb (Cat. No. RMK0269)
Standard	31.25-2000 pg/mL	Recombinant Mouse VCAM-1/CD106 Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.