

Rabbit anti-Mouse TPO mAb

Catalog No.: RMK0308

Basic Information

Catalog No.
RMK0308

Category
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
/

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Iodination and coupling of the hormonogenic tyrosines in thyroglobulin to yield the thyroid hormones T3 and T4

Immunogen Information

Immunogen
Recombinant Mouse TPO

Cross-Reactivity


/

Assay Applications

Mouse TPO Sandwich Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4 µg/mL	Rabbit anti-Mouse TPO mAb (Cat. No. RMK0308)
ELISA Detection	0.0075-0.03 µg/mL	Rabbit anti-Mouse TPO mAb (Cat. No. RMK0309)
Standard	31.25-2000 pg/mL	Recombinant Mouse TPO Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.