

Anti- Mouse TPO mAb (CAP)

Catalog No.: RMK0308

Basic Information

Catalog No.
RMK0308

Category
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS with 0.05%ProClin 300,PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Iodination and coupling of the hormonogenic tyrosines in thyroglobulin to yield the thyroid hormones T3 and T4

Immunogen Information

Immunogen
Recombinant Mouse TPO

Cross-Reactivity


/

Assay Applications

Mouse TPO Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4 ug/mL	Rabbit anti-Mouse TPO (CAP)(Cat. No.RMK0307)
ELISA Detection	0.0075-0.03ug/mL	Rabbit anti-Mouse TPO (DET)(Cat. No.RMK0308)
Standard	31.25-2000pg/mL	Recombinant Mouse TPO Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.