

# Anti- Mouse TPO mAb (DET)

**Catalog No.:** RMK0309

## Basic Information

**Catalog No.**  
RMK0309

**Category**  
Elisa Antibody Kit

**Application**  
ELISA

## Product Information

**Ig Type**  
Rabbit IgG


**Purification**  
Affinity purification


**Endotoxin Level**

**Storage**  
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2um filtered solution in PBS with 0.05%ProClin 300,PH 7.4.

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

Iodination and coupling of the hormonogenic tyrosines in thyroglobulin to yield the thyroid hormones T3 and T4

## Immunogen Information

**Immunogen**  
Recombinant Mouse TPO

## Cross-Reactivity

/


## Assay Applications

Mouse TPO Sandwich ELISA Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4 ug/mL	Rabbit anti-Mouse TPO (CAP)(Cat. No.RMK0307)
ELISA Detection	0.0075-0.03ug/mL	Rabbit anti-Mouse TPO (DET)(Cat. No.RMK0308)
Standard	31.25-2000pg/mL	Recombinant Mouse TPO Protein

## Validation Data

---



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.