

Anti-Mouse TPO mAb (DET)

Catalog No.: RMK0309

Basic Information

Catalog No.

RMK0309

Catagory

Elisa Antibody Kit

Application

ELISA

Product Information

Ig Type

Rabbit IgG

Purification

Affinity purification

Endotoxin Level

/

Storage


Store at -20°C.


Avoid repeated freeze-thaw cycles.

Formulation

Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This gene encodes a humoral growth factor necessary for megakaryocyte proliferation and maturation, as well as for thrombopoiesis. The encoded protein is a ligand for the product of the myeloproliferative leukemia virus oncogene. Alternative splicing results in multiple transcript variants.

Immunogen Information

Immunogen

Recombinant Mouse TPO

Cross-Reactivity


/

Assay Applications

Mouse TPO Sandwich Immunoassay

	Recommended Concentration	Sample
ELISA Capture	1-4 µg/mL	Rabbit anti-Mouse TPO mAb (Cat. No. RMK0308)
ELISA Detection	0.0075-0.03 µg/mL	Rabbit anti-Mouse TPO mAb (Cat. No. RMK0309)
Standard	31.25-2000 pg/mL	Recombinant Mouse TPO Protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.