

Rabbit anti-Mouse GM-CSF mAb(DET)

Catalog No.: RMK0322

Basic Information

Catalog No.
RMK0322

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Mouse GM-CSF Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

GM-CSF was initially characterized as a factor that can support the in vitro colony formation of granulocyte-macrophage progenitors. It is also a growth factor for erythroid, megakaryocyte, and eosinophil progenitors. GM-CSF is produced by a number of different cell types (including T cells, B cells, macrophages, mast cells, endothelial cells, fibroblasts, and adipocytes) in response to cytokine or inflammatory stimuli

Immunogen Information

Immunogen
Recombinant Mouse GM-CSF Protein

Cross-Reactivity

/

Assay Applications

Mouse GM-CSF multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Mouse GM-CSF mAb(Cat. No.RMK0321)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Mouse GM-CSF mAb(Cat. No.RMK0322)
Standard	6.86-5000pg/mL	Recombinant-Mouse GM-CSF Protein

Validation Data

Standard Curve Chart

