

Rabbit anti-Mouse CXCL13 mAb(CAP)

Catalog No.: RMK0328

Basic Information

Catalog No.
RMK0328

Category
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Mouse CXCL13 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

CXCL13, also known as B-lymphocyte chemoattractant (BLC), is a CXC chemokine that is constitutively expressed in secondary lymphoid organs. Mouse BCA-1 cDNA encodes a precursor protein of 109 amino acid residues with a putative leader sequence of 21 residues. Mature mouse BCA-1 shares 64% amino acid sequence similarity with the human protein and 23-34% amino acid sequence identity with other known CXC chemokines

Immunogen Information

Immunogen
Recombinant Mouse CXCL13 Protein

Cross-Reactivity

/

Assay Applications

Mouse CXCL13 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Mouse CXCL13 mAb(Cat. No.RMK0328)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Mouse CXCL13 mAb(Cat. No.RMK0329)
Standard	6.86-5000pg/mL	Recombinant-Mouse CXCL13 Protein

Validation Data

Standard Curve Chart

