

Rabbit anti-Mouse IL-13 mAb (DET)

Catalog No.: RMK0343

Basic Information

Catalog No.
RMK0343

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative[]0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Mouse IL-13 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

IL-13 (Interleukin-13) plays a key role in the pathogenesis of allergy, cancer, and tissue fibrosis. It is secreted by several helper T cell subsets, NK cells, mast cells, eosinophils, basophils, and visceral smooth muscle cells. IL-13 suppresses the production of proinflammatory cytokines and other cytotoxic substances by macrophages, fibroblasts, and endothelial cells.

Immunogen Information

Immunogen
Recombinant Mouse IL-13 Protein

Cross-Reactivity

/

Assay Applications

Mouse IL-13 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Mouse IL-13 mAb(Cat. No.RMK0342)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Mouse IL-13 mAb(Cat. No.RMK0343)
Standard	6.86-5000pg/mL	Recombinant-Mouse IL-13 Protein

Validation Data

Standard Curve Chart

