

Rabbit anti-Human IgM mAb(DET)

Catalog No.: RMK0354

Basic Information

Catalog No.
RMK0354

Category
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human IgM Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Our NorthernLights fluorescent secondary antibodies are bright and resistant to photobleaching. We are currently offering secondary antibodies recognizing mouse, rat, goat, sheep, and rabbit IgG as well as chicken IgY. These reagents are available with three distinct excitation and emission maxima, making them ideal for multi-color fluorescence microscopy.

Immunogen Information

Immunogen
Human IgM

Cross-Reactivity

/

Assay Applications

Human IgM multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human IgM mAb(Cat. No.RMK0353)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human IgM mAb(Cat. No.RMK0354)
Standard	6.86-5000pg/mL	Recombinant Human IgM Protein

Validation Data

Standard Curve Chart

