

Rabbit anti- Human TSLP mAb(DET)

Catalog No.: RMK0357

Basic Information

Catalog No.
RMK0357

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative[]0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human TSLP Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Thymic Stromal Lymphopoietin (TSLP) was originally identified as an activity from the conditioned medium of a mouse thymic stromal cell line that promoted the development of B cells (1 - 3). The activities of mouse TSLP overlap with, but are distinct from, those of mouse IL-7. Both mouse TSLP and IL-7 can co-stimulate growth of thymocytes and mature T cells, and support B lymphopoiesis in long-term cultures of fetal liver cells and bone-marrow cells

Immunogen Information

Immunogen
Recombinant Human TSLP Protein

Cross-Reactivity

/

Assay Applications

Human TSLP multiplex assay

| | Recommended Concentration | Sample |
|---------------------|---------------------------|---|
| Multiplex Capture | 3-20ug/mL | Rabbit anti-Human TSLP mAb(Cat. No.RMK0356) |
| Multiplex Detection | 0.017-2ug/mL | Rabbit anti-Human TSLP mAb(Cat. No.RMK0357) |
| Standard | 6.86-5000pg/mL | Recombinant Human TSLP Protein |

Validation Data

Standard Curve Chart

