

# Rabbit anti-Human TNF RII/TNFRSF1B mAb(DET)

**Catalog No.:** RMK0366

## Basic Information

**Catalog No.**  
RMK0366

**Category**  
Elisa Antibody Kit

**Application**  
multiplex assay

## Product Information

**Ig Type**  
Rabbit IgG

**Purification**  
Affinity purification


**Endotoxin Level**

**Storage**  
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human TNF RII/TNFRSF1B Antibody

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

TNF RII, also known as p75 or p80 TNF R, is a prototypic member of the TNF receptor superfamily and has been designated TNFRSF1B. It binds TNF with high affinity and is expressed primarily in cells of the immune system.

## Immunogen Information

**Immunogen**  
Recombinant Human TNF RII/TNFRSF1B Protein

## Cross-Reactivity

/

## Assay Applications

Human TNF RII/TNFRSF1B multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human TNF RII/TNFRSF1B mAb(Cat. No.RMK0365)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human TNF RII/TNFRSF1B mAb(Cat. No.RMK0366)
Standard	6.86-5000pg/mL	Recombinant Human TNF RII/TNFRSF1B Protein

## Validation Data

---

Standard Curve Chart

