

# Rabbit anti-Human TFF3 mAb(CAP)

**Catalog No.:** RMK0368

## Basic Information

**Catalog No.**  
RMK0368

**Catagory**  
Elisa Antibody Kit

**Application**  
multiplex assay

## Product Information

**Ig Type**  
Rabbit IgG


**Purification**  
Affinity purification


**Endotoxin Level**

**Storage**  
Store at -20°C. Avoid freeze / thaw cycles. Preservative[]0.05%ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human TFF3 Antibody

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

TFF3 (Trefoil factor 3; also hP1.B and intestinal trefoil factor) is a 7 kDa member of the trefoil factor family of protease-resistant peptides. It is expressed by multiple epithelial cell types, acts in an autocrine or paracrine manner and contributes to epithelial cell migration and survival

## Immunogen Information

**Immunogen**  
Recombinant Human TFF3 Protein

## Cross-Reactivity

/

## Assay Applications

Human TFF3 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human TFF3 mAb(Cat. No.RMK0368)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human TFF3 mAb(Cat. No.RMK0369)
Standard	6.86-5000pg/mL	Recombinant Human TFF3 Protein

## Validation Data

---

Standard Curve Chart

