

Rabbit anti-Human E-Selectin/CD62E mAb(CAP)

Catalog No.: RMK0370

Basic Information

Catalog No.
RMK0370

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative[]0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human E-Selectin/CD62E Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

E-Selectin, also known as Endothelial Leukocyte Adhesion Molecule 1 (ELAM-1), and Leukocyte-Endothelial Cell Adhesion Molecule 2 (LECAM-2), has been designated CD62E. It is a cell-surface glycoprotein expressed by activated endothelium and it mediates the adhesion of blood neutrophils.

Immunogen Information

Immunogen
Recombinant Human E-Selectin/CD62E Protein

Cross-Reactivity

/

Assay Applications

Human E-Selectin/CD62E multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human TE-Selectin/CD62E mAb(Cat. No.RMK0370)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human E-Selectin/CD62E mAb(Cat. No.RMK0371)
Standard	6.86-5000pg/mL	Recombinant Human E-Selectin/CD62E Protein

Validation Data

Standard Curve Chart

