

# Rabbit anti-Human E-Selectin/CD62E mAb(DET)

**Catalog No.:** RMK0371

## Basic Information

**Catalog No.**  
RMK0371

**Catagory**  
Elisa Antibody Kit

**Application**  
multiplex assay

## Product Information

**Ig Type**  
Rabbit IgG

**Purification**  
Affinity purification

**Endotoxin Level**

**Storage**  
Store at -20°C. Avoid freeze / thaw cycles. Preservative[]0.05%ProClin 300.  
**Avoid repeated freeze-thaw cycles.**

**Formulation**  
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human E-Selectin/CD62E Antibody

## Contact

 | [order@abclonal.com](mailto:order@abclonal.com)

 | [support@abclonal.com](mailto:support@abclonal.com)

 | [www.abclonal.com](http://www.abclonal.com)

## Background

E-Selectin, also known as Endothelial Leukocyte Adhesion Molecule 1 (ELAM-1), and Leukocyte-Endothelial Cell Adhesion Molecule 2 (LECAM-2), has been designated CD62E. It is a cell-surface glycoprotein expressed by activated endothelium and it mediates the adhesion of blood neutrophils.

## Immunogen Information

**Immunogen**  
Recombinant Human E-Selectin/CD62E Protein

## Cross-Reactivity

/

## Assay Applications

Human E-Selectin/CD62E multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human TE-Selectin/CD62E mAb(Cat. No.RMK0370)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human E-Selectin/CD62E mAb(Cat. No.RMK0371)
Standard	6.86-5000pg/mL	Recombinant Human E-Selectin/CD62E Protein

# Validation Data

---

Standard Curve Chart

