

Rabbit anti-Human ICAM-1/CD54 mAb(CAP)

Catalog No.: RMK0383

Basic Information

Catalog No.
RMK0383

Category
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human ICAM-1/CD54 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

ICAM-1 is a member of the immunoglobulin superfamily whose expression is upregulated on leukocytes, epithelial cells and resting endothelial cells in response to inflammatory signals. ICAM-1 binds the leukocyte integrins LFA-1 and Mac-1.

Immunogen Information

Immunogen
Recombinant Human ICAM-1/CD54 Protein

Cross-Reactivity

/

Assay Applications

Human ICAM-1/CD54 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human ICAM-1/CD54 mAb(Cat. No.RMK0383)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human ICAM-1/CD54 mAb(Cat. No.RMK0384)
Standard	6.86-5000pg/mL	Recombinant Human ICAM-1/CD54 Protein

Validation Data

Standard Curve Chart

