

Rabbit anti-Human ICAM-1/CD54 mAb

Catalog No.: RMK0383

Basic Information

Catalog No.
RMK0383

Category
Elisa Antibody Kit

Application
Multiplex

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
N/A

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.

Immunogen Information

Immunogen
Recombinant Human ICAM-1/CD54 Protein

Cross-Reactivity

/

Assay Applications

Human ICAM-1/CD54 Sandwich Immunoassay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human ICAM-1/CD54 mAb(Cat. No.RMK0383)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human ICAM-1/CD54 mAb(Cat. No.RMK0384)
Standard	6.86-5000pg/mL	Recombinant Human ICAM-1/CD54 Protein

Validation Data

Standard Curve Chart

