

Rabbit anti-Human CXCL9 mAb(DET)

Catalog No.: RMK0389

Basic Information

Catalog No.
RMK0389

Category
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human CXCL9 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

CXCL9, also known as MIG, is a chemokine that is induced in macrophages and in CNS glial cells in response to IFN-gamma. CXCL9 signals through CXCR3 to induce the chemoattraction of activated T cells and tumor-infiltrating lymphocytes.

Immunogen Information

Immunogen
Recombinant Human CXCL9 Protein

Cross-Reactivity

/

Assay Applications

Human CXCL9 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human CXCL9 mAb(Cat. No.RMK0388)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human CXCL9 mAb(Cat. No.RMK0389)
Standard	6.86-5000pg/mL	Recombinant Human CXCL9 Protein

Validation Data

Standard Curve Chart

