

Rabbit anti-Human LRG1 mAb(CAP)

Catalog No.: RMK0390

Basic Information

Catalog No.
RMK0390

Category
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human LRG1 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

The leucine-rich repeat (LRR) family of proteins, including LRG1, have been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. LRG1 is expressed during granulocyte differentiation (O'Donnell et al., 2002 [PubMed 12223515]).

Immunogen Information

Immunogen
Recombinant Human LRG1 Protein

Cross-Reactivity

/

Assay Applications

Human LRG1 multiplex assay

| | Recommended Concentration | Sample |
|---------------------|---------------------------|---|
| Multiplex Capture | 3-20ug/mL | Rabbit anti-Human LRG1 mAb(Cat. No.RMK0390) |
| Multiplex Detection | 0.017-2ug/mL | Rabbit anti-Human LRG1 mAb(Cat. No.RMK0291) |
| Standard | 6.86-5000pg/mL | Recombinant Human LRG1 Protein |

Validation Data

Standard Curve Chart

