

Rabbit anti-Human LRG1 mAb(CAP)

Catalog No.: RMK0390

Basic Information

Catalog No.
RMK0390

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative[]0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human LRG1 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

LRG-1 (Leucine-Rich alpha-2-Glycoprotein 1; also LRHG [in rodent]) is a 45-55 kDa, monomeric, variably glycosylated, secreted protein member of the leucine-rich repeat family of molecules. It should not be confused with either LRIG-1/ Leucine-Rich repeats and Ig-like domains-1, or LRG1/LIM-RhoGap protein 1

Immunogen Information

Immunogen
Recombinant Human LRG1 Protein

Cross-Reactivity

/

Assay Applications

Human LRG1 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human LRG1 mAb(Cat. No.RMK0390)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human LRG1 mAb(Cat. No.RMK0291)
Standard	6.86-5000pg/mL	Recombinant Human LRG1 Protein

Validation Data

Standard Curve Chart

