

Rabbit anti-Human TRAIL-R3/CD263 mAb (DET)

Catalog No.: RMK0404

Basic Information

Catalog No.
RMK0404

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative[]0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human TRAIL-R3/CD263 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor contains an extracellular TRAIL-binding domain and a transmembrane domain, but no cytoplasmic death domain. This receptor is not capable of inducing apoptosis, and is thought to function as an antagonistic receptor that protects cells from TRAIL-induced apoptosis. This gene was found to be a p53-regulated DNA damage-inducible gene. The expression of this gene was detected in many normal tissues but not in most cancer cell lines, which may explain the specific sensitivity of cancer cells to the apoptosis-inducing activity of TRAIL.

Immunogen Information

Immunogen
Recombinant Human TRAIL-R3/CD263 Protein

Cross-Reactivity

/

Assay Applications

Human TRAIL-R3/CD263 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human TRAIL-R3/CD263 mAb(Cat. No.RMK0403)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human TRAIL-R3/CD263 mAb(Cat. No.RMK0404)
Standard	6.86-5000pg/mL	Recombinant Human TRAIL-R3/CD263 Protein

Validation Data

Standard Curve Chart

