

Rabbit anti-Human TRAIL-R3/CD263 mAb (DET)

Catalog No.: RMK0404

Basic Information

Catalog No.
RMK0404

Category
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05% ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS, PH 7.4 containing Human TRAIL-R3/CD263 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Human TRAIL R3, also called DcR1 (decoy receptor 1), LIT, and TRID, is a glycosyl-phosphatidylinositol-linked membrane protein which binds TRAIL (Apo2 Ligand) with high affinity. In the new TNF superfamily nomenclature, TRAIL R3 is referred to as TNFRSF10C. TRAIL R3 has the TRAIL-binding extracellular cysteine-rich domains but lacks the intracellular signalling domain. As a result, binding of TRAIL to TRAIL R3 does not transduce an apoptosis signal

Immunogen Information

Immunogen
Recombinant Human TRAIL-R3/CD263 Protein

Cross-Reactivity

/

Assay Applications

Human TRAIL-R3/CD263 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human TRAIL-R3/CD263 mAb(Cat. No.RMK0403)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human TRAIL-R3/CD263 mAb(Cat. No.RMK0404)
Standard	6.86-5000pg/mL	Recombinant Human TRAIL-R3/CD263 Protein

Validation Data

Standard Curve Chart

