

Rabbit anti-Human Peroxiredoxin 5 (PRDX5) mAb (CAP)

Catalog No.: RMK0405

Basic Information

Catalog No.
RMK0405

Category
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05% ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS, PH 7.4 containing Human Peroxiredoxin 5 Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Peroxiredoxin 5 (Prx-5; also known as AOEB166 and thioredoxin reductase) is a 22 kDa, widely expressed mitochondrial antioxidant enzyme that belongs to the atypical 2-Cys class of the TSA/ahpC family of peroxiredoxins

Immunogen Information

Immunogen
Recombinant Human Peroxiredoxin 5 Protein

Cross-Reactivity

/

Assay Applications

Human PRDX5 multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human PRDX5 mAb(Cat. No.RMK0405)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human PRDX5 mAb(Cat. No.RM18052)
Standard	6.86-5000pg/mL	Recombinant Human PRDX5 Protein

Validation Data

Standard Curve Chart

