

Rabbit anti-Human sGSN mAb (CAP)

Catalog No.: RMK0408

Basic Information

Catalog No.
RMK0408

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG


Purification
Affinity purification


Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human sGSN Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Gelsolin, also known as GSN, actin-depolymerizing factor/ADF, AGEL and Brevin, is a 90-95 kDa member of the villin/gelsolin family. Widely expressed, Gelsolin binds to actin and fibronectin, and is found both secreted in plasma and in cytoplasm. Cytoplasmic Gelsolin lacks 51 N-terminal amino acids (aa) present in the secreted protein

Immunogen Information

Immunogen
Recombinant Human sGSN Protein

Cross-Reactivity

/

Assay Applications

Human sGSN multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human sGSN mAb(Cat. No.RMK0408)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human sGSN mAb(Cat. No.RMK0409)
Standard	6.86-5000pg/mL	Recombinant Human sGSN Protein

Validation Data

Standard Curve Chart

