

Rabbit anti-Human sGSN mAb (DET)

Catalog No.: RMK0409

Basic Information

Catalog No.
RMK0409

Catagory
Elisa Antibody Kit

Application
multiplex assay

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C. Avoid freeze / thaw cycles. Preservative 0.05%ProClin 300.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2um filtered solution in PBS,PH 7.4 containing Human sGSN Antibody

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

The protein encoded by this gene binds to the "plus" ends of actin monomers and filaments to prevent monomer exchange. The encoded calcium-regulated protein functions in both assembly and disassembly of actin filaments. Defects in this gene are a cause of familial amyloidosis Finnish type (FAF). Multiple transcript variants encoding several different isoforms have been found for this gene.

Immunogen Information

Immunogen
Recombinant Human sGSN Protein

Cross-Reactivity

/

Assay Applications

Human sGSN multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human sGSN mAb(Cat. No.RMK0408)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human sGSN mAb(Cat. No.RMK0409)
Standard	6.86-5000pg/mL	Recombinant Human sGSN Protein

Validation Data

Standard Curve Chart

