

Rabbit anti-Human CXCL10/IP-10 mAb

Catalog No.: RMK0528

Basic Information

Catalog No.
RMK0528

Catagory
Elisa Antibody Kit

Application
Multiplex

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2µm filtered solution in PBS with 0.05% Proclin300, PH 7.4.

Contact

		order@abclonal.com
		support@abclonal.com
		www.abclonal.com

Background

This antimicrobial gene encodes a chemokine of the CXC subfamily and ligand for the receptor CXCR3. Binding of this protein to CXCR3 results in pleiotropic effects, including stimulation of monocytes, natural killer and T-cell migration, and modulation of adhesion molecule expression. This gene may also be a key regulator of the 'cytokine storm' immune response to SARS-CoV-2 infection.

Immunogen Information

Immunogen
Recombinant Human CXCL10/IP-10 Protein

Cross-Reactivity

/

Assay Applications

Human CXCL10/IP-10 Sandwich Immunoassay

		Recommended Concentration	Sample
Multiplex	Capture	1.2-4.8ug/mL	Rabbit anti-Human CXCL10/IP-10 mAb(Cat. No. RMK0527)
	Detection	0.125-0.5ug/mL	Rabbit anti-Human CXCL10/IP-10 mAb(Cat. No. RMK0528)
	Standard	6.86-5000pg/mL	Recombinant Human CXCL10/IP-10 Protein(Cat. No. RP01638)