

Mouse anti-Human IgA mAb

Catalog No.: RMK0574

Basic Information

Catalog No.
RMK0574

Catagory
Elisa Antibody Kit

Application
ELISA

Product Information

Ig Type
Mouse IgG2b

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, PH 7.4.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Contributes to immunoglobulin receptor binding activity. Involved in antibacterial humoral response; glomerular filtration; and positive regulation of respiratory burst. Located in extracellular space. Part of monomeric IgA immunoglobulin complex and secretory dimeric IgA immunoglobulin complex.

Immunogen Information

Immunogen

Cross-Reactivity

/

Assay Applications