

Rabbit anti-Human Granzyme B mAb

Catalog No.: RMK0585

Basic Information

Catalog No.
RMK0585

Category
Elisa Antibody Kit

Application
Multiplex

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

This gene encodes a member of the granzyme subfamily of proteins, part of the peptidase S1 family of serine proteases. The encoded preproprotein is secreted by natural killer (NK) cells and cytotoxic T lymphocytes (CTLs) and proteolytically processed to generate the active protease, which induces target cell apoptosis. This protein also processes cytokines and degrades extracellular matrix proteins, and these roles are implicated in chronic inflammation and wound healing. Expression of this gene may be elevated in human patients with cardiac fibrosis.

Immunogen Information

Immunogen
Recombinant Human Granzyme B protein

Cross-Reactivity


Human

Assay Applications

Human Granzyme B Sandwich Immunoassay

		Recommended Concentration	Sample
Multiplex	Capture	3-20µg/mL	Rabbit anti-Human Granzyme B mAb(Cat. No. RMK0584)
	Detection	0.017-2 µg/mL	Rabbit anti-Human Granzyme B mAb(Cat. No. RMK0585)
	Standard	6.86-5000 pg/mL	Recombinant Human Granzyme B protein

Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.