

Rabbit anti-Mouse CXCL10/IP-10 mAb

Catalog No.: RMK0606

Basic Information

Catalog No.
RMK0606

Category
Elisa Antibody Kit

Application

Product Information

Ig Type
Rabbit IgG

Purification
Affinity purification

Endotoxin Level
N/A

Storage
Store at -20°C.
Avoid repeated freeze-thaw cycles.

Formulation
Supplied as a 0.2 µm filtered solution in PBS with 0.05% Proclin 300, pH 7.3.

Contact

 | order@abclonal.com

 | support@abclonal.com

 | www.abclonal.com

Background

Predicted to enable chemoattractant activity; chemokine receptor binding activity; and heparin binding activity. Involved in cellular response to interleukin-17. Acts upstream of or within several processes, including defense response to virus; negative regulation of myoblast differentiation; and negative regulation of myoblast fusion. Located in external side of plasma membrane and extracellular space. Is expressed in several structures, including alimentary system; axial skeleton; hemolymphoid system; lacrimal gland; and liver. Human ortholog(s) of this gene implicated in hepatitis B; middle cerebral artery infarction; and type 1 diabetes mellitus. Orthologous to human CXCL10 (C-X-C motif chemokine ligand 10).

Immunogen Information

Immunogen
Recombinant Mouse CXCL10 protein

Cross-Reactivity

Mouse

Assay Applications
